

CE 1434




Science&Care


# Novacheck® SARS-CoV-2 Antigen Rapid Test Selftest



Assay Format	Lateral flow test Immuno-chromatographic in vitro test
Test Type	Qualitative
Instrument	Device-free
Sample Material	Nasal
Target Antigen	Nucleocapsid (N)
Readout time	15 minutes
Sensitivity	95.91%
Specificity	99.9%
Storage temperature	2- 30 °C



### Screening




SARS- CoV-2, also known as the Covid-19 virus, causes an acute respiratory infectious disease. The main source of infection is currently people infected with the virus, even those who have an asymptomatic course. According to recent epidemiological studies, the incubation period varies from 1 to 14 days, but mostly from 3 to 7 days.

Pandemic can only be contained by interrupting the transmission of infection. Considering that most transmission is caused by asymptomatic virus carriers, systematic and repeated screening of the population has a very special role. Without screening, we deprive ourselves of the opportunity to prevent further de-escalation.



### Novacheck® SARS-CoV-2 Antigen Rapid Test

The Novacheck® SARS-CoV-2 Antigen Rapid Test is an immunochromatographic assay for the qualitative detection of SARS-CoV-2 nucleocapsid antigen in human nasal swabs. For nasal sampling, the sample is collected from the anterior region of the nose.




### Clinical Validation

The SARS-COV-2 rapid test kit was performed on 342 nasal swabs collected from symptomatic patients who had symptom onset within 7 days. A limited number of patients with symptom onset for more than 7 days and asymptomatic patients were included in the clinical study (n=342). The sample size was relatively significant, positive agreement 95.91% (328/342) and negative agreement 99.9% (517/517). The test is designed for lay user.

Method	PCR			
	Result	positive	negative	Final result
Novacheck	positive	328	0	328
	negative	14	517	531
<b>Final result</b>		<b>342</b>	<b>517</b>	<b>859</b>
Sensitivity:	95.91% [328/342]	(95% CI: 93.86-99.04)		
Specificity:	99.9% [517/517]	(95% CI: 98.53-100.00)		
Accuracy:	98.37% (517+328) + (328+0+14+517)			

### How does Lateral Flow Antigen Rapid Test work?



### Eurosurveillance Evaluation



Experts from the Paul Ehrlich Institute, in collaboration with researchers from other institutions, examined a total of 122 COVID-19 antigen rapid tests for sensitivity and thus for their ability to detect the SARS-CoV-2 virus.



**The result:**  
The quality of the tests varied widely. The Novacheck® SARS-CoV-2 Antigen Rapid Test was among 96 rapid antigen tests that met the required criteria; 26 tests did not provide the required sensitivity.

  
**CERTIFICATE**

**EC Certificate No. 1434-IVDD-009/2022**  
EC Design-examination  
Directive 98/79/EC concerning  
in vitro diagnostic medical devices

Polish Centre for Testing and Certification certifies  
that manufactured by:

**NOVATECH TIBBI CIHAZ ÜRÜNLERİ SANAYİ VE TİCARET A.Ş.**  
2.Organize Sanayi Bölgesi Hacı Sani Konukoğlu Bulvarı  
83228 Nolu, Sokak No:17 Sehitkamil / Gaziantep, Turkey

in vitro diagnostic medical devices  
for self-testing

**Novacheck® SARS-CoV-2 Antigen Rapid Test MY-28**

In terms of design documentation, comply with requirements  
of Annex I8 (Section 8) to Directive 98/79/EC (as amended)  
implemented into Polish law,  
as evidenced by the audit conducted by the PCBC

Validity of the Certificate: from 23.01.2022 to 27.05.2025  
The date of expiry of the certificate: 27.05.2022  
The date of the first issue of the Certificate: 23.01.2022

**CE 1434**

Issued under the Contract No. AN-156/2021  
Application No. 2723/2021  
Certificate issue: the audit program:  
Version: 22.01.2022  
Module A1

Digitally signed  
by Aleksandra  
Kostrzewa  
President

POLISH CENTRE FOR TESTING AND CERTIFICATION | 02 664 89 000, 400 P.Ławiecka Street, 01-483 48 41 286, e-mail: info@pcbc.gov.pl

**Technical Universal Verification**  
**CERTIFICATE**

This Certificate has been issued to:  
KOMPLEKSYWNE LABORATORIUM DIAGNOSTYCZNE T.C.M.T. A.S.

ISO 13485 : 2016  
For the scope of activities described in the certificate

CSCD - HASIŁO PRZEKAZANE W TRYBIE INFORMACYJNEJ

NAC

**Technical Universal Verification**  
**CERTIFICATE**

This Certificate has been issued to:  
KOMPLEKSYWNE LABORATORIUM DIAGNOSTYCZNE T.C.M.T. A.S.

ISO 9001 : 2015  
For the scope of activities described in the certificate

CSCD - HASIŁO PRZEKAZANE W TRYBIE INFORMACYJNEJ

NAC

**Technical Universal Verification**  
**CERTIFICATE**

This Certificate has been issued to:  
KOMPLEKSYWNE LABORATORIUM DIAGNOSTYCZNE T.C.M.T. A.S.

ISO 22716 - GMP  
For the scope of activities described in the certificate

CSCD - HASIŁO PRZEKAZANE W TRYBIE INFORMACYJNEJ

NAC



**PEI**  
evaluiert

EU common list



Science&Care

Science & Care  
Vertriebs GmbH

Science & Care Vertriebs GmbH  
Südring 59, D-69514 Laudendach  
+49 6201 99910  
info@sciencecare.de  
www.sciencecare.de