

Sample ID	Specimen type for the antigen test	Specimen collection date and time for the antigen test		Specimen testing date and time for the antigen test		Antigen test result with the analyzer or reader value	Specimen type for Reference KIT	Specimen collection date and time for Reference KIT		Specimen testing date and time for Reference KIT		Name of Reference KIT used as the comparator	RT-PCR test results (+/-)	RT-PCR TEST VALUE RESULT (Ct VALUES)	Number of days post-onset of patient symptoms	Patient age and gender, if available
		Date	Time	Date	Time			Date	Time	Date	Time					
SMP-217	Nasal	4/18/20	9:17	4/18/20	9:17	+	Nasopharynx	4/18/20	9:17	4/18/20	15:00	Bio-tech	+	18,56	1	M
SMP-859	Nasal	30/08/2020	10:53	30/08/2020	10:53	+	Nasopharynx	30/08/2020	10:53	30/08/2020	15:00	Bio-tech	+	19,01	1	FM
SMP-124	Nasal	5/04/2020	09:18	5/04/2020	09:18	+	Nasopharynx	5/04/2020	09:18	5/04/2020	15:00	Bio-tech	+	19,80	1	FM
SMP-627	Nasal	27/07/2020	09:17	27/07/2020	09:17	+	Nasopharynx	27/07/2020	09:17	27/07/2020	15:00	Bio-tech	+	20,16	1	FM
SMP-502	Nasal	7/06/2020	09:49	7/06/2020	09:49	+	Nasopharynx	7/06/2020	09:49	7/06/2020	15:00	Bio-tech	+	20,34	2	M
SMP-509	Nasal	14/06/2020	09:14	14/06/2020	09:14	+	Nasopharynx	14/06/2020	09:14	14/06/2020	15:00	Bio-tech	+	20,34	2	M
SMP-562	Nasal	12/07/2020	10:00	12/07/2020	10:00	+	Nasopharynx	12/07/2020	10:00	12/07/2020	15:00	Bio-tech	+	20,34	1	M
SMP-799	Nasal	24/08/2020	09:54	24/08/2020	09:54	+	Nasopharynx	24/08/2020	09:54	24/08/2020	15:00	Bio-tech	+	20,38	1	M
SMP-614	Nasal	26/07/2020	09:18	26/07/2020	09:18	+	Nasopharynx	26/07/2020	09:18	26/07/2020	15:00	Bio-tech	+	20,43	1	FM
SMP-448	Nasal	5/25/20	11:40	5/25/20	11:40	+	Nasopharynx	5/25/20	11:40	5/25/20	15:00	Bio-tech	+	20,45	2	M
SMP-741	Nasal	19/08/2020	11:32	19/08/2020	11:32	+	Nasopharynx	19/08/2020	11:32	19/08/2020	15:00	Bio-tech	+	20,46	2	M
SMP-240	Nasal	4/21/20	10:02	4/21/20	10:02	+	Nasopharynx	4/21/20	10:02	4/21/20	15:00	Bio-tech	+	20,54	1	M
SMP-561	Nasal	11/07/2020	09:41	11/07/2020	09:41	+	Nasopharynx	11/07/2020	09:41	11/07/2020	15:00	Bio-tech	+	20,65	1	M
SMP-810	Nasal	24/08/2020	10:03	24/08/2020	10:03	+	Nasopharynx	24/08/2020	10:03	24/08/2020	15:00	Bio-tech	+	20,65	1	FM
SMP-660	Nasal	2/08/2020	10:02	2/08/2020	10:02	+	Nasopharynx	2/08/2020	10:02	2/08/2020	15:00	Bio-tech	+	20,96	1	M
SMP-295	Nasal	4/29/20	09:48	4/29/20	09:48	+	Nasopharynx	4/29/20	09:48	4/29/20	15:00	Bio-tech	+	20,99	2	FM
SMP-439	Nasal	5/23/20	09:14	5/23/20	09:14	+	Nasopharynx	5/23/20	09:14	5/23/20	15:00	Bio-tech	+	21,01	2	M
SMP-153	Nasal	8/04/2020	11:05	8/04/2020	11:05	+	Nasopharynx	8/04/2020	11:05	8/04/2020	15:00	Bio-tech	+	21,08	2	FM
SMP-271	Nasal	4/25/20	09:49	4/25/20	09:49	+	Nasopharynx	4/25/20	09:49	4/25/20	15:00	Bio-tech	+	21,08	3	FM
SMP-752	Nasal	21/08/2020	9:33	21/08/2020	9:33	+	Nasopharynx	21/08/2020	9:33	21/08/2020	15:00	Bio-tech	+	21,13	2	FM
SMP-470	Nasal	5/28/20	9:44	5/28/20	9:44	+	Nasopharynx	5/28/20	9:44	5/28/20	15:00	Bio-tech	+	21,21	3	M
SMP-197	Nasal	4/15/20	11:11	4/15/20	11:11	+	Nasopharynx	4/15/20	11:11	4/15/20	15:00	Bio-tech	+	21,34	3	M
SMP-491	Nasal	5/30/20	9:56	5/30/20	9:56	+	Nasopharynx	5/30/20	9:56	5/30/20	15:00	Bio-tech	+	21,34	1	M
SMP-589	Nasal	20/07/2020	09:54	20/07/2020	09:54	+	Nasopharynx	20/07/2020	09:54	20/07/2020	15:00	Bio-tech	+	21,55	1	M
SMP-185	Nasal	4/13/20	10:12	4/13/20	10:12	+	Nasopharynx	4/13/20	10:12	4/13/20	15:00	Bio-tech	+	21,58	3	M
SMP-156	Nasal	8/04/2020	11:48	8/04/2020	11:48	+	Nasopharynx	8/04/2020	11:48	8/04/2020	15:00	Bio-tech	+	21,72	2	M
SMP-429	Nasal	5/21/20	10:33	5/21/20	10:33	+	Nasopharynx	5/21/20	10:33	5/21/20	15:00	Bio-tech	+	21,72	2	FM
SMP-831	Nasal	26/08/2020	11:12	26/08/2020	11:12	+	Nasopharynx	26/08/2020	11:12	26/08/2020	15:00	Bio-tech	+	21,78	1	M
SMP-190	Nasal	4/14/20	09:24	4/14/20	09:24	+	Nasopharynx	4/14/20	09:24	4/14/20	15:00	Bio-tech	+	21,87	3	M
SMP-210	Nasal	4/17/20	09:35	4/17/20	09:35	+	Nasopharynx	4/17/20	09:35	4/17/20	15:00	Bio-tech	+	21,87	3	M
SMP-159	Nasal	9/04/2020	09:14	9/04/2020	09:14	+	Nasopharynx	9/04/2020	09:14	9/04/2020	15:00	Bio-tech	+	21,91	2	M
SMP-179	Nasal	12/04/2020	11:36	12/04/2020	11:36	+	Nasopharynx	12/04/2020	11:36	12/04/2020	15:00	Bio-tech	+	21,91	3	M
SMP-119	Nasal	4/04/2020	09:30	4/04/2020	09:30	+	Nasopharynx	4/04/2020	09:30	4/04/2020	15:00	Bio-tech	+	22,09	2	FM
SMP-79	Nasal	25/03/2020	10:33	25/03/2020	10:33	+	Nasopharynx	25/03/2020	10:33	25/03/2020	15:00	Bio-tech	+	22,11	2	FM
SMP-334	Nasal	4/05/2020	11:18	4/05/2020	11:18	+	Nasopharynx	4/05/2020	11:18	4/05/2020	15:00	Bio-tech	+	22,11	1	FM
SMP-161	Nasal	9/04/2020	10:00	9/04/2020	10:00	+	Nasopharynx	9/04/2020	10:00	9/04/2020	15:00	Bio-tech	+	22,16	3	M
SMP-01	Nasal	21/03/2020	8:46	21/03/2020	8:46	+	Nasopharynx	21/03/2020	8:46	21/03/2020	15:00	Bio-tech	+	22,20	1	M
SMP-343	Nasal	6/05/2020	8:04	6/05/2020	8:04	+	Nasopharynx	6/05/2020	8:04	6/05/2020	15:00	Bio-tech	+	22,27	1	M
SMP-184	Nasal	4/13/20	9:43	4/13/20	9:43	+	Nasopharynx	4/13/20	9:43	4/13/20	15:00	Bio-tech	+	22,32	2	M
SMP-431	Nasal	5/21/20	11:27	5/21/20	11:27	+	Nasopharynx	5/21/20	11:27	5/21/20	15:00	Bio-tech	+	22,32	3	FM
SMP-595	Nasal	22/07/2020	08:47	22/07/2020	08:47	+	Nasopharynx	22/07/2020	08:47	22/07/2020	15:00	Bio-tech	+	22,33	4	FM
SMP-811	Nasal	24/08/2020	10:32	24/08/2020	10:32	+	Nasopharynx	24/08/2020	10:32	24/08/2020	15:00	Bio-tech	+	22,33	2	M
SMP-512	Nasal	16/06/2020	10:04	16/06/2020	10:04	+	Nasopharynx	16/06/2020	10:04	16/06/2020	15:00	Bio-tech	+	22,34	2	FM
SMP-445	Nasal	5/25/20	09:55	5/25/20	09:55	+	Nasopharynx	5/25/20	09:55	5/25/20	15:00	Bio-tech	+	22,39	3	M
SMP-192	Nasal	4/14/20	10:33	4/14/20	10:33	+	Nasopharynx	4/14/20	10:33	4/14/20	15:00	Bio-tech	+	22,42	2	FM
SMP-275	Nasal	4/26/20	10:53	4/26/20	10:53	+	Nasopharynx	4/26/20	10:53	4/26/20	15:00	Bio-tech	+	22,42	2	M
SMP-771	Nasal	22/08/2020	09:41	22/08/2020	09:41	+	Nasopharynx	22/08/2020	09:41	22/08/2020	15:00	Bio-tech	+	22,43	1	M
SMP-457	Nasal	5/26/20	09:44	5/26/20	09:44	+	Nasopharynx	5/26/20	09:44	5/26/20	15:00	Bio-tech	+	22,46	2	M
SMP-276	Nasal	4/26/20	11:10	4/26/20	11:10	+	Nasopharynx	4/26/20	11:10	4/26/20	15:00	Bio-tech	+	22,50	3	FM
SMP-840	Nasal	28/08/2020	09:44	28/08/2020	09:44	+	Nasopharynx	28/08/2020	09:44	28/08/2020	15:00	Bio-tech	+	22,54	1	M
SMP-651	Nasal	1/08/2020	10:00	1/08/2020	10:00	+	Nasopharynx	1/08/2020	10:00	1/08/2020	15:00	Bio-tech	+	22,55	3	M
SMP-527	Nasal	29/06/2020	09:44	29/06/2020	09:44	+	Nasopharynx	29/06/2020	09:44	29/06/2020	15:00	Bio-tech	+	22,56	2	M
SMP-148	Nasal	8/04/2020	9:21	8/04/2020	9:21	+	Nasopharynx	8/04/2020	9:21	8/04/2020	15:00	Bio-tech	+	22,57	3	M
SMP-166	Nasal	10/04/2020	10:33	10/04/2020	10:33	+	Nasopharynx	10/04/2020	10:33	10/04/2020	15:00	Bio-tech	+	22,57	3	M
SMP-531	Nasal	1/07/2020	10:32	1/07/2020	10:32	+	Nasopharynx	1/07/2020	10:32	1/07/2020	15:00	Bio-tech	+	22,58	2	M
SMP-483	Nasal	5/29/20	10:32	5/29/20	10:32	+	Nasopharynx	5/29/20	10:32	5/29/20	15:00	Bio-tech	+	22,65	1	M
SMP-209	Nasal	4/16/20	11:52	4/16/20	11:52	+	Nasopharynx	4/16/20	11:52	4/16/20	15:00	Bio-tech	+	22,68	2	FM
SMP-534	Nasal	2/07/2020	10:43	2/07/2020	10:43	+	Nasopharynx	2/07/2020	10:43	2/07/2020	15:00	Bio-tech	+	22,68	2	M
SMP-557	Nasal	10/07/2020	09:17	10/07/2020	09:17	+	Nasopharynx	10/07/2020	09:17	10/07/2020	15:00	Bio-tech	+	22,68	1	FM
SMP-200	Nasal	4/15/20	11:43	4/15/20	11:43	+	Nasopharynx	4/15/20	11:43	4/15/20	15:00	Bio-tech	+	22,72	2	FM
SMP-279	Nasal	4/27/20	09:23	4/27/20	09:23	+	Nasopharynx	4/27/20	09:23	4/27/20	15:00	Bio-tech	+	22,73	2	FM
SMP-497	Nasal	6/3/20	09:17	6/3/20	09:17	+	Nasopharynx	6/3/20	09:17	6/3/20	15:00	Bio-tech	+	22,82	1	M
SMP-657	Nasal	2/08/2020	9:56	2/08/2020	9:56	+	Nasopharynx	2/08/2020	9:56	2/08/2020	15:00	Bio-tech	+	22,82	3	M
SMP-672	Nasal	6/08/2020	10:25	6/08/2020	10:25	+	Nasopharynx	6/08/2020	10:25	6/08/2020	15:00	Bio-tech	+	22,83	3	M
SMP-167	Nasal	11/04/2020	09:33	11/04/2020	09:33	+	Nasopharynx	11/04/2020	09:33	11/04/2020	15:00	Bio-tech	+	22,85	3	M
SMP-539	Nasal	5/07/2020	09:30	5/07/2020	09:30	+	Nasopharynx	5/07/2020	09:30	5/07/2020	15:00	Bio-tech	+	22,99	2	M
SMP-12	Nasal	23/03/2020	09:22	23/03/2020	09:22	+	Nasopharynx	23/03/2020	09:22	23/03/2020	15:00	Bio-tech	+	23,01	2	FM
SMP-157	Nasal	9/04/2020	08:47	9/04/2020	08:47	+	Nasopharynx	9/04/2020	08:47	9/04/2020	15:00	Bio-tech	+	23,01	3	M
SMP-402	Nasal	5/15/20	10:33	5/15/20	10:33	+	Nasopharynx	5/15/20	10:33	5/15/20	15:00	Bio-tech	+	23,01	2	FM
SMP-680	Nasal	7/08/2020	09:24	7/08/2020	09:24	+	Nasopharynx	7/08/2020	09:24	7/08/2020	15:00	Bio-tech	+	23,04	3	M

SMP-207	Nasal	4/16/20	11:16	4/16/20	11:16	+	Nasopharynx	4/16/20	11:16	4/16/20	15:00	Bio-tech	+	23,05	1	FM
SMP-133	Nasal	6/04/2020	10:23	6/04/2020	10:23	+	Nasopharynx	6/04/2020	10:23	6/04/2020	15:00	Bio-tech	+	23,10	1	M
SMP-37	Nasal	14/03/2020	09:44	14/03/2020	09:44	+	Nasopharynx	14/03/2020	09:44	14/03/2020	15:00	Bio-tech	+	23,11	2	M
SMP-407	Nasal	5/17/20	09:11	5/17/20	09:11	+	Nasopharynx	5/17/20	09:11	5/17/20	15:00	Bio-tech	+	23,11	3	M
SMP-781	Nasal	23/08/2020	09:27	23/08/2020	09:27	+	Nasopharynx	23/08/2020	09:27	23/08/2020	15:00	Bio-tech	+	23,11	3	M
SMP-786	Nasal	23/08/2020	11:05	23/08/2020	11:05	+	Nasopharynx	23/08/2020	11:05	23/08/2020	15:00	Bio-tech	+	23,11	2	M
SMP-793	Nasal	24/08/2020	10:27	24/08/2020	10:27	+	Nasopharynx	24/08/2020	10:27	24/08/2020	15:00	Bio-tech	+	23,11	2	M
SMP-801	Nasal	24/08/2020	9:13	24/08/2020	9:13	+	Nasopharynx	24/08/2020	9:13	24/08/2020	15:00	Bio-tech	+	23,11	3	M
SMP-804	Nasal	24/08/2020	09:03	24/08/2020	09:03	+	Nasopharynx	24/08/2020	09:03	24/08/2020	15:00	Bio-tech	+	23,11	2	M
SMP-724	Nasal	17/08/2020	09:13	17/08/2020	09:13	+	Nasopharynx	17/08/2020	09:13	17/08/2020	15:00	Bio-tech	+	23,28	2	M
SMP-751	Nasal	21/08/2020	8:58	21/08/2020	8:58	+	Nasopharynx	21/08/2020	8:58	21/08/2020	15:00	Bio-tech	+	23,28	4	M
SMP-624	Nasal	27/07/2020	9:46	27/07/2020	9:46	+	Nasopharynx	27/07/2020	9:46	27/07/2020	15:00	Bio-tech	+	23,32	4	FM
SMP-440	Nasal	5/23/20	09:38	5/23/20	09:38	+	Nasopharynx	5/23/20	09:38	5/23/20	15:00	Bio-tech	+	23,33	3	M
SMP-438	Nasal	5/23/20	09:02	5/23/20	09:02	+	Nasopharynx	5/23/20	09:02	5/23/20	15:00	Bio-tech	+	23,38	2	M
SMP-348	Nasal	6/05/2020	09:49	6/05/2020	09:49	+	Nasopharynx	6/05/2020	09:49	6/05/2020	15:00	Bio-tech	+	23,44	4	M
SMP-486	Nasal	5/30/20	09:10	5/30/20	09:10	+	Nasopharynx	5/30/20	09:10	5/30/20	15:00	Bio-tech	+	23,45	3	FM
SMP-121	Nasal	5/04/2020	09:58	5/04/2020	09:58	+	Nasopharynx	5/04/2020	09:58	5/04/2020	15:00	Bio-tech	+	23,46	4	M
SMP-219	Nasal	4/18/20	10:33	4/18/20	10:33	+	Nasopharynx	4/18/20	10:33	4/18/20	15:00	Bio-tech	+	23,53	2	FM
SMP-180	Nasal	12/04/2020	11:47	12/04/2020	11:47	+	Nasopharynx	12/04/2020	11:47	12/04/2020	15:00	Bio-tech	+	23,54	1	FM
SMP-469	Nasal	5/28/20	09:30	5/28/20	09:30	+	Nasopharynx	5/28/20	09:30	5/28/20	15:00	Bio-tech	+	23,55	2	FM
SMP-852	Nasal	29/08/2020	09:43	29/08/2020	09:43	+	Nasopharynx	29/08/2020	09:43	29/08/2020	15:00	Bio-tech	+	23,55	2	FM
SMP-691	Nasal	10/08/2020	09:49	10/08/2020	09:49	+	Nasopharynx	10/08/2020	09:49	10/08/2020	15:00	Bio-tech	+	23,62	3	FM
SMP-222	Nasal	4/18/20	11:49	4/18/20	11:49	+	Nasopharynx	4/18/20	11:49	4/18/20	15:00	Bio-tech	+	23,65	2	M
SMP-241	Nasal	4/21/20	10:21	4/21/20	10:21	+	Nasopharynx	4/21/20	10:21	4/21/20	15:00	Bio-tech	+	23,65	3	M
SMP-715	Nasal	15/08/2020	09:48	15/08/2020	09:48	+	Nasopharynx	15/08/2020	09:48	15/08/2020	15:00	Bio-tech	+	23,65	2	FM
SMP-829	Nasal	26/08/2020	10:32	26/08/2020	10:32	+	Nasopharynx	26/08/2020	10:32	26/08/2020	15:00	Bio-tech	+	23,66	2	M
SMP-22	Nasal	26/03/2020	10:04	26/03/2020	10:04	+	Nasopharynx	26/03/2020	10:04	26/03/2020	15:00	Bio-tech	+	23,71	2	FM
SMP-375	Nasal	10/05/2020	8:12	10/05/2020	8:12	+	Nasopharynx	10/05/2020	8:12	10/05/2020	15:00	Bio-tech	+	23,77	3	M
SMP-779	Nasal	23/08/2020	10:33	23/08/2020	10:33	+	Nasopharynx	23/08/2020	10:33	23/08/2020	15:00	Bio-tech	+	23,77	2	FM
SMP-301	Nasal	4/29/20	13:12	4/29/20	13:12	+	Nasopharynx	4/29/20	13:12	4/29/20	15:00	Bio-tech	+	23,87	3	M
SMP-772	Nasal	22/08/2020	10:00	22/08/2020	10:00	+	Nasopharynx	22/08/2020	10:00	22/08/2020	15:00	Bio-tech	+	23,87	1	M
SMP-338	Nasal	5/05/2020	8:56	5/05/2020	8:56	+	Nasopharynx	5/05/2020	8:56	5/05/2020	15:00	Bio-tech	+	23,88	2	M
SMP-204	Nasal	4/16/20	9:46	4/16/20	9:46	+	Nasopharynx	4/16/20	9:46	4/16/20	15:00	Bio-tech	+	23,89	4	FM
SMP-225	Nasal	4/19/20	09:48	4/19/20	09:48	+	Nasopharynx	4/19/20	09:48	4/19/20	15:00	Bio-tech	+	23,92	2	FM
SMP-229	Nasal	4/20/20	9:17	4/20/20	9:17	+	Nasopharynx	4/20/20	9:17	4/20/20	15:00	Bio-tech	+	23,96	2	M
SMP-410	Nasal	5/17/20	09:43	5/17/20	09:43	+	Nasopharynx	5/17/20	09:43	5/17/20	15:00	Bio-tech	+	24,02	2	FM
SMP-250	Nasal	4/22/20	11:47	4/22/20	11:47	+	Nasopharynx	4/22/20	11:47	4/22/20	15:00	Bio-tech	+	24,03	1	FM
SMP-548	Nasal	8/07/2020	09:17	8/07/2020	09:17	+	Nasopharynx	8/07/2020	09:17	8/07/2020	15:00	Bio-tech	+	24,03	2	M
SMP-453	Nasal	5/25/20	12:55	5/25/20	12:55	+	Nasopharynx	5/25/20	12:55	5/25/20	15:00	Bio-tech	+	24,04	4	M
SMP-309	Nasal	30/04/2020	12:14	30/04/2020	12:14	+	Nasopharynx	30/04/2020	12:14	30/04/2020	15:00	Bio-tech	+	24,05	1	FM
SMP-232	Nasal	4/20/20	11:04	4/20/20	11:04	+	Nasopharynx	4/20/20	11:04	4/20/20	15:00	Bio-tech	+	24,07	2	M
SMP-373	Nasal	9/05/2020	10:27	9/05/2020	10:27	+	Nasopharynx	9/05/2020	10:27	9/05/2020	15:00	Bio-tech	+	24,07	2	M
SMP-418	Nasal	5/19/20	11:44	5/19/20	11:44	+	Nasopharynx	5/19/20	11:44	5/19/20	15:00	Bio-tech	+	24,12	4	M
SMP-359	Nasal	7/05/2020	11:33	7/05/2020	11:33	+	Nasopharynx	7/05/2020	11:33	7/05/2020	15:00	Bio-tech	+	24,12	2	FM
SMP-92	Nasal	1/04/2020	10:04	1/04/2020	10:04	+	Nasopharynx	1/04/2020	10:04	1/04/2020	15:00	Bio-tech	+	24,13	2	FM
SMP-388	Nasal	12/05/2020	10:24	12/05/2020	10:24	+	Nasopharynx	12/05/2020	10:24	12/05/2020	15:00	Bio-tech	+	24,17	2	M
SMP-759	Nasal	21/08/2020	11:45	21/08/2020	11:45	+	Nasopharynx	21/08/2020	11:45	21/08/2020	15:00	Bio-tech	+	24,22	4	FM
SMP-318	Nasal	1/05/2020	12:24	1/05/2020	12:24	+	Nasopharynx	1/05/2020	12:24	1/05/2020	15:00	Bio-tech	+	24,24	2	M
SMP-555	Nasal	10/07/2020	11:02	10/07/2020	11:02	+	Nasopharynx	10/07/2020	11:02	10/07/2020	15:00	Bio-tech	+	24,26	2	M
SMP-287	Nasal	4/28/20	09:17	4/28/20	09:17	+	Nasopharynx	4/28/20	09:17	4/28/20	15:00	Bio-tech	+	24,31	1	M
SMP-471	Nasal	5/28/20	09:58	5/28/20	09:58	+	Nasopharynx	5/28/20	09:58	5/28/20	15:00	Bio-tech	+	24,31	4	M
SMP-244	Nasal	4/21/20	11:03	4/21/20	11:03	+	Nasopharynx	4/21/20	11:03	4/21/20	15:00	Bio-tech	+	24,32	2	M
SMP-611	Nasal	25/07/2020	09:58	25/07/2020	09:58	+	Nasopharynx	25/07/2020	09:58	25/07/2020	15:00	Bio-tech	+	24,32	4	M
SMP-247	Nasal	4/22/20	09:44	4/22/20	09:44	+	Nasopharynx	4/22/20	09:44	4/22/20	15:00	Bio-tech	+	24,40	2	M
SMP-115	Nasal	4/04/2020	11:12	4/04/2020	11:12	+	Nasopharynx	4/04/2020	11:12	4/04/2020	15:00	Bio-tech	+	24,41	3	M
SMP-709	Nasal	14/08/2020	10:33	14/08/2020	10:33	+	Nasopharynx	14/08/2020	10:33	14/08/2020	15:00	Bio-tech	+	24,41	2	FM
SMP-173	Nasal	11/04/2020	10:55	11/04/2020	10:55	+	Nasopharynx	11/04/2020	10:55	11/04/2020	15:00	Bio-tech	+	24,44	4	FM
SMP-367	Nasal	8/05/2020	11:47	8/05/2020	11:47	+	Nasopharynx	8/05/2020	11:47	8/05/2020	15:00	Bio-tech	+	24,44	3	M
SMP-778	Nasal	23/08/2020	10:21	23/08/2020	10:21	+	Nasopharynx	23/08/2020	10:21	23/08/2020	15:00	Bio-tech	+	24,44	3	M
SMP-699	Nasal	12/08/2020	09:23	12/08/2020	09:23	+	Nasopharynx	12/08/2020	09:23	12/08/2020	15:00	Bio-tech	+	24,52	2	FM
SMP-41	Nasal	15/03/2020	10:32	15/03/2020	10:32	+	Nasopharynx	15/03/2020	10:32	15/03/2020	15:00	Bio-tech	+	24,55	2	M
SMP-139	Nasal	7/04/2020	9:44	7/04/2020	9:44	+	Nasopharynx	7/04/2020	9:44	7/04/2020	15:00	Bio-tech	+	24,55	2	FM
SMP-270	Nasal	4/25/20	09:43	4/25/20	09:43	+	Nasopharynx	4/25/20	09:43	4/25/20	15:00	Bio-tech	+	24,55	2	FM
SMP-400	Nasal	5/15/20	09:24	5/15/20	09:24	+	Nasopharynx	5/15/20	09:24	5/15/20	15:00	Bio-tech	+	24,55	3	M
SMP-800	Nasal	24/08/2020	10:02	24/08/2020	10:02	+	Nasopharynx	24/08/2020	10:02	24/08/2020	15:00	Bio-tech	+	24,55	1	M
SMP-300	Nasal	4/29/20	12:38	4/29/20	12:38	+	Nasopharynx	4/29/20	12:38	4/29/20	15:00	Bio-tech	+	24,56	3	M
SMP-695	Nasal	10/08/2020	10:53	10/08/2020	10:53	+	Nasopharynx	10/08/2020	10:53	10/08/2020	15:00	Bio-tech	+	24,60	2	M
SMP-685	Nasal	8/08/2020	08:57	8/08/2020	08:57	+	Nasopharynx	8/08/2020	08:57	8/08/2020	15:00	Bio-tech	+	24,64	2	M
SMP-254	Nasal	4/23/20	10:43	4/23/20	10:43	+	Nasopharynx	4/23/20	10:43	4/23/20	15:00	Bio-tech	+	24,66	2	M
SMP-344	Nasal	6/05/2020	8:32	6/05/2020	8:32	+	Nasopharynx	6/05/2020	8:32	6/05/2020	15:00	Bio-tech	+	24,67	4	FM
SMP-447	Nasal	5/25/20	11:03	5/25/20	11:03	+	Nasopharynx	5/25/20	11:03	5/25/20	15:00	Bio-tech	+	24,67	3	M
SMP-633	Nasal	28/07/2020	10:11	28/07/2020	10:11	+	Nasopharynx	28/07/2020	10:11	28/07/2020	15:00	Bio-tech	+	24,67	4	M
SMP-671	Nasal	6/08/2020	10:32	6/08/2020	10:32	+	Nasopharynx	6/08/2020	10:32	6/08/2020	15:00	Bio-tech	+	24,72	2	M
SMP-257	Nasal	4/23/20	11:37	4/23/20	11:37	+	Nasopharynx	4/23/20	11:37	4/23/20	15:00	Bio-tech	+	24,78	3	FM
SMP-415	Nasal	5/18/20	10:53	5/18/20	10:53	+	Nasopharynx	5/18/20	10:53	5/18/20	15:00	Bio-tech	+	24,78	2	M
SMP-421	Nasal	5/20/20	09:41	5/20/20	09:41	+	Nasopharynx	5/20/20	09:41							

SMP-645	Nasal	31/07/2020	09:48	31/07/2020	09:48	+	Nasopharynx	31/07/2020	09:48	31/07/2020	15:00	Bio-tech	+	24,78	2	M
SMP-650	Nasal	1/08/2020	10:38	1/08/2020	10:38	+	Nasopharynx	1/08/2020	10:38	1/08/2020	15:00	Bio-tech	+	24,78	3	M
SMP-795	Nasal	24/08/2020	11:12	24/08/2020	11:12	+	Nasopharynx	24/08/2020	11:12	24/08/2020	15:00	Bio-tech	+	24,78	3	M
SMP-236	Nasal	4/21/20	8:33	4/21/20	8:33	+	Nasopharynx	4/21/20	8:33	4/21/20	15:00	Bio-tech	+	24,84	3	M
SMP-494	Nasal	5/30/20	10:35	5/30/20	10:35	+	Nasopharynx	5/30/20	10:35	5/30/20	15:00	Bio-tech	+	24,87	4	M
SMP-783	Nasal	23/08/2020	10:38	23/08/2020	10:38	+	Nasopharynx	23/08/2020	10:38	23/08/2020	15:00	Bio-tech	+	24,88	2	M
SMP-60	Nasal	21/03/2020	09:43	21/03/2020	09:43	+	Nasopharynx	21/03/2020	09:43	21/03/2020	15:00	Bio-tech	+	24,89	2	M
SMP-141	Nasal	7/04/2020	10:41	7/04/2020	10:41	+	Nasopharynx	7/04/2020	10:41	7/04/2020	15:00	Bio-tech	+	24,89	1	M
SMP-432	Nasal	5/22/20	09:49	5/22/20	09:49	+	Nasopharynx	5/22/20	09:49	5/22/20	15:00	Bio-tech	+	24,89	2	M
SMP-24	Nasal	27/03/2020	09:13	27/03/2020	09:13	+	Nasopharynx	27/03/2020	09:13	27/03/2020	15:00	Bio-tech	+	24,91	2	M
SMP-422	Nasal	5/20/20	10:00	5/20/20	10:00	+	Nasopharynx	5/20/20	10:00	5/20/20	15:00	Bio-tech	+	24,91	1	M
SMP-670	Nasal	6/08/2020	10:03	6/08/2020	10:03	+	Nasopharynx	6/08/2020	10:03	6/08/2020	15:00	Bio-tech	+	24,91	1	FM
SMP-576	Nasal	16/07/2020	11:05	16/07/2020	11:05	+	Nasopharynx	16/07/2020	11:05	16/07/2020	15:00	Bio-tech	+	24,92	2	M
SMP-319	Nasal	1/05/2020	13:12	1/05/2020	13:12	-	Nasopharynx	1/05/2020	13:12	1/05/2020	15:00	Bio-tech	+	24,12	3	M
False <25Ct						1										
Total correct						159							160			
SMP-23	Nasal	26/03/2020	10:40	26/03/2020	10:40	+	Nasopharynx	26/03/2020	10:40	26/03/2020	15:00	Bio-tech	+	25,06	2	M
SMP-322	Nasal	2/05/2020	10:45	2/05/2020	10:45	+	Nasopharynx	2/05/2020	10:45	2/05/2020	15:00	Bio-tech	+	25,07	3	M
SMP-291	Nasal	4/28/20	11:27	4/28/20	11:27	+	Nasopharynx	4/28/20	11:27	4/28/20	15:00	Bio-tech	+	25,09	3	M
SMP-520	Nasal	23/06/2020	10:02	23/06/2020	10:02	+	Nasopharynx	23/06/2020	10:02	23/06/2020	15:00	Bio-tech	+	25,10	1	M
SMP-61	Nasal	21/03/2020	09:49	21/03/2020	09:49	+	Nasopharynx	21/03/2020	09:49	21/03/2020	15:00	Bio-tech	+	25,12	3	FM
SMP-572	Nasal	15/07/2020	09:49	15/07/2020	09:49	+	Nasopharynx	15/07/2020	09:49	15/07/2020	15:00	Bio-tech	+	25,12	2	FM
SMP-601	Nasal	24/07/2020	10:32	24/07/2020	10:32	+	Nasopharynx	24/07/2020	10:32	24/07/2020	15:00	Bio-tech	+	25,12	2	M
SMP-643	Nasal	30/07/2020	10:38	30/07/2020	10:38	+	Nasopharynx	30/07/2020	10:38	30/07/2020	15:00	Bio-tech	+	25,13	2	M
SMP-259	Nasal	4/24/20	09:30	4/24/20	09:30	+	Nasopharynx	4/24/20	09:30	4/24/20	15:00	Bio-tech	+	25,18	2	FM
SMP-728	Nasal	18/08/2020	9:32	18/08/2020	9:32	+	Nasopharynx	18/08/2020	9:32	18/08/2020	15:00	Bio-tech	+	25,18	2	M
SMP-639	Nasal	29/07/2020	10:33	29/07/2020	10:33	+	Nasopharynx	29/07/2020	10:33	29/07/2020	15:00	Bio-tech	+	25,23	2	FM
SMP-578	Nasal	16/07/2020	09:02	16/07/2020	09:02	+	Nasopharynx	16/07/2020	09:02	16/07/2020	15:00	Bio-tech	+	25,25	2	M
SMP-702	Nasal	13/08/2020	9:00	13/08/2020	9:00	+	Nasopharynx	13/08/2020	9:00	13/08/2020	15:00	Bio-tech	+	25,25	1	M
SMP-341	Nasal	5/05/2020	09:49	5/05/2020	09:49	+	Nasopharynx	5/05/2020	09:49	5/05/2020	15:00	Bio-tech	+	25,27	3	FM
SMP-629	Nasal	27/07/2020	09:23	27/07/2020	09:23	+	Nasopharynx	27/07/2020	09:23	27/07/2020	15:00	Bio-tech	+	25,30	2	FM
SMP-622	Nasal	27/07/2020	9:23	27/07/2020	9:23	+	Nasopharynx	27/07/2020	9:23	27/07/2020	15:00	Bio-tech	+	25,31	2	M
SMP-91	Nasal	1/04/2020	10:00	1/04/2020	10:00	+	Nasopharynx	1/04/2020	10:00	1/04/2020	15:00	Bio-tech	+	25,32	3	M
SMP-836	Nasal	27/08/2020	10:14	27/08/2020	10:14	+	Nasopharynx	27/08/2020	10:14	27/08/2020	15:00	Bio-tech	+	25,33	3	M
SMP-846	Nasal	28/08/2020	13:22	28/08/2020	13:22	+	Nasopharynx	28/08/2020	13:22	28/08/2020	15:00	Bio-tech	+	25,33	4	FM
SMP-317	Nasal	1/05/2020	11:44	1/05/2020	11:44	+	Nasopharynx	1/05/2020	11:44	1/05/2020	15:00	Bio-tech	+	25,34	2	M
SMP-592	Nasal	21/07/2020	09:37	21/07/2020	09:37	+	Nasopharynx	21/07/2020	09:37	21/07/2020	15:00	Bio-tech	+	25,34	3	FM
SMP-304	Nasal	30/04/2020	11:13	30/04/2020	11:13	+	Nasopharynx	30/04/2020	11:13	30/04/2020	15:00	Bio-tech	+	25,36	2	M
SMP-737	Nasal	19/08/2020	09:44	19/08/2020	09:44	+	Nasopharynx	19/08/2020	09:44	19/08/2020	15:00	Bio-tech	+	25,36	2	M
SMP-615	Nasal	26/07/2020	08:57	26/07/2020	08:57	+	Nasopharynx	26/07/2020	08:57	26/07/2020	15:00	Bio-tech	+	25,43	2	M
SMP-762	Nasal	22/08/2020	10:23	22/08/2020	10:23	+	Nasopharynx	22/08/2020	10:23	22/08/2020	15:00	Bio-tech	+	25,44	2	M
SMP-770	Nasal	22/08/2020	09:35	22/08/2020	09:35	+	Nasopharynx	22/08/2020	09:35	22/08/2020	15:00	Bio-tech	+	25,45	3	M
SMP-280	Nasal	4/27/20	09:35	4/27/20	09:35	+	Nasopharynx	4/27/20	09:35	4/27/20	15:00	Bio-tech	+	25,48	3	M
SMP-818	Nasal	25/08/2020	11:36	25/08/2020	11:36	+	Nasopharynx	25/08/2020	11:36	25/08/2020	15:00	Bio-tech	+	25,51	5	M
SMP-612	Nasal	25/07/2020	10:33	25/07/2020	10:33	+	Nasopharynx	25/07/2020	10:33	25/07/2020	15:00	Bio-tech	+	25,52	2	FM
SMP-82	Nasal	30/03/2020	09:49	30/03/2020	09:49	+	Nasopharynx	30/03/2020	09:49	30/03/2020	15:00	Bio-tech	+	25,55	2	M
SMP-94	Nasal	1/04/2020	09:13	1/04/2020	09:13	+	Nasopharynx	1/04/2020	09:13	1/04/2020	15:00	Bio-tech	+	25,55	2	M
SMP-127	Nasal	6/04/2020	8:35	6/04/2020	8:35	+	Nasopharynx	6/04/2020	8:35	6/04/2020	15:00	Bio-tech	+	25,55	3	M
SMP-196	Nasal	4/15/20	10:04	4/15/20	10:04	+	Nasopharynx	4/15/20	10:04	4/15/20	15:00	Bio-tech	+	25,55	4	M
SMP-339	Nasal	5/05/2020	09:25	5/05/2020	09:25	+	Nasopharynx	5/05/2020	09:25	5/05/2020	15:00	Bio-tech	+	25,55	4	FM
SMP-355	Nasal	7/05/2020	10:50	7/05/2020	10:50	+	Nasopharynx	7/05/2020	10:50	7/05/2020	15:00	Bio-tech	+	25,55	4	M
SMP-370	Nasal	8/05/2020	12:38	8/05/2020	12:38	+	Nasopharynx	8/05/2020	12:38	8/05/2020	15:00	Bio-tech	+	25,55	3	M
SMP-411	Nasal	5/17/20	09:49	5/17/20	09:49	+	Nasopharynx	5/17/20	09:49	5/17/20	15:00	Bio-tech	+	25,55	3	M
SMP-515	Nasal	19/06/2020	09:55	19/06/2020	09:55	+	Nasopharynx	19/06/2020	09:55	19/06/2020	15:00	Bio-tech	+	25,55	3	M
SMP-522	Nasal	24/06/2020	09:37	24/06/2020	09:37	+	Nasopharynx	24/06/2020	09:37	24/06/2020	15:00	Bio-tech	+	25,55	3	FM
SMP-523	Nasal	26/06/2020	09:55	26/06/2020	09:55	+	Nasopharynx	26/06/2020	09:55	26/06/2020	15:00	Bio-tech	+	25,55	4	M
SMP-237	Nasal	4/21/20	8:56	4/21/20	8:56	+	Nasopharynx	4/21/20	8:56	4/21/20	15:00	Bio-tech	+	25,56	3	M
SMP-551	Nasal	10/07/2020	09:49	10/07/2020	09:49	+	Nasopharynx	10/07/2020	09:49	10/07/2020	15:00	Bio-tech	+	25,58	3	FM
SMP-126	Nasal	5/04/2020	09:04	5/04/2020	09:04	+	Nasopharynx	5/04/2020	09:04	5/04/2020	15:00	Bio-tech	+	25,59	4	M
SMP-616	Nasal	26/07/2020	10:04	26/07/2020	10:04	+	Nasopharynx	26/07/2020	10:04	26/07/2020	15:00	Bio-tech	+	25,65	4	M
SMP-791	Nasal	24/08/2020	10:00	24/08/2020	10:00	+	Nasopharynx	24/08/2020	10:00	24/08/2020	15:00	Bio-tech	+	25,65	3	M
SMP-20	Nasal	26/03/2020	09:38	26/03/2020	09:38	+	Nasopharynx	26/03/2020	09:38	26/03/2020	15:00	Bio-tech	+	25,66	3	M
SMP-792	Nasal	24/08/2020	10:04	24/08/2020	10:04	+	Nasopharynx	24/08/2020	10:04	24/08/2020	15:00	Bio-tech	+	25,66	2	FM
SMP-814	Nasal	25/08/2020	10:43	25/08/2020	10:43	+	Nasopharynx	25/08/2020	10:43	25/08/2020	15:00	Bio-tech	+	25,66	2	M
SMP-214	Nasal	4/17/20	10:43	4/17/20	10:43	+	Nasopharynx	4/17/20	10:43	4/17/20	15:00	Bio-tech	+	25,67	4	M
SMP-598	Nasal	23/07/2020	10:24	23/07/2020	10:24	+	Nasopharynx	23/07/2020	10:24	23/07/2020	15:00	Bio-tech	+	25,67	2	M
SMP-726	Nasal	17/08/2020	10:33	17/08/2020	10:33	+	Nasopharynx	17/08/2020	10:33	17/08/2020	15:00	Bio-tech	+	25,67	3	M
SMP-62	Nasal	21/03/2020	10:23	21/03/2020	10:23	+	Nasopharynx	21/03/2020	10:23	21/03/2020	15:00	Bio-tech	+	25,70	2	M
SMP-544	Nasal	6/07/2020	09:18	6/07/2020	09:18	+	Nasopharynx	6/07/2020	09:18	6/07/2020	15:00	Bio-tech	+	25,70	1	FM
SMP-740	Nasal	19/08/2020	11:03	19/08/2020	11:03	+	Nasopharynx	19/08/2020	11:03	19/08/2020	15:00	Bio-tech	+	25,70	1	FM
SMP-735	Nasal	19/08/2020	08:47	19/08/2020	08:47	+	Nasopharynx	19/08/2020	08:47	19/08/2020	15:00	Bio-tech	+	25,72	4	FM
SMP-584	Nasal	18/07/2020	10:13	18/07/2020	10:13	+	Nasopharynx	18/07/2020	10:13	18/07/2020	15:00	Bio-tech	+	25,77	2	M
SMP-480	Nasal	5/29/20	09:43	5/29/20	09:43	+	Nasopharynx	5/29/20	09:43	5/29/20	15:00	Bio-tech	+	25,79	2	FM
SMP-437	Nasal	5/22/20	12:47	5/22/20	12:47	+	Nasopharynx	5/22/20	12:47	5/22/20	15:00	Bio-tech	+	25,87	3	M
SMP-632	Nasal	28/07/2020	10:00	28/07/2020	10:00	+	Nasopharynx	28/07/2020	10:00	28/07/2020	15:00	Bio-tech	+	26,01	1	M
SMP-583	Nasal	18/07/2020	9:40	18/07/2020	9:40	+	Nasopharynx									

SMP-501	Nasal	6/06/2020	09:27	6/06/2020	09:27	+	Nasopharynx	6/06/2020	09:27	6/06/2020	15:00	Bio-tech	+	26,10	3	M
SMP-570	Nasal	15/07/2020	09:22	15/07/2020	09:22	+	Nasopharynx	15/07/2020	09:22	15/07/2020	15:00	Bio-tech	+	26,11	5	M
SMP-842	Nasal	28/08/2020	11:24	28/08/2020	11:24	+	Nasopharynx	28/08/2020	11:24	28/08/2020	15:00	Bio-tech	+	26,11	3	M
SMP-25	Nasal	27/03/2020	09:55	27/03/2020	09:55	+	Nasopharynx	27/03/2020	09:55	27/03/2020	15:00	Bio-tech	+	26,19	3	M
SMP-462	Nasal	5/26/20	12:25	5/26/20	12:25	+	Nasopharynx	5/26/20	12:25	5/26/20	15:00	Bio-tech	+	26,20	3	M
SMP-780	Nasal	23/08/2020	09:22	23/08/2020	09:22	+	Nasopharynx	23/08/2020	09:22	23/08/2020	15:00	Bio-tech	+	26,21	5	M
SMP-822	Nasal	25/08/2020	13:32	25/08/2020	13:32	+	Nasopharynx	25/08/2020	13:32	25/08/2020	15:00	Bio-tech	+	26,23	3	M
SMP-296	Nasal	4/29/20	11:05	4/29/20	11:05	+	Nasopharynx	4/29/20	11:05	4/29/20	15:00	Bio-tech	+	26,26	2	M
SMP-559	Nasal	10/07/2020	09:23	10/07/2020	09:23	+	Nasopharynx	10/07/2020	09:23	10/07/2020	15:00	Bio-tech	+	26,27	2	FM
SMP-507	Nasal	12/06/2020	08:47	12/06/2020	08:47	+	Nasopharynx	12/06/2020	08:47	12/06/2020	15:00	Bio-tech	+	26,30	3	M
SMP-52	Nasal	19/03/2020	10:33	19/03/2020	10:33	+	Nasopharynx	19/03/2020	10:33	19/03/2020	15:00	Bio-tech	+	26,31	2	FM
SMP-140	Nasal	7/04/2020	10:00	7/04/2020	10:00	+	Nasopharynx	7/04/2020	10:00	7/04/2020	15:00	Bio-tech	+	26,31	3	M
SMP-81	Nasal	25/03/2020	09:27	25/03/2020	09:27	+	Nasopharynx	25/03/2020	09:27	25/03/2020	15:00	Bio-tech	+	26,32	3	M
SMP-608	Nasal	25/07/2020	09:44	25/07/2020	09:44	+	Nasopharynx	25/07/2020	09:44	25/07/2020	15:00	Bio-tech	+	26,33	5	M
SMP-129	Nasal	6/04/2020	09:25	6/04/2020	09:25	+	Nasopharynx	6/04/2020	09:25	6/04/2020	15:00	Bio-tech	+	26,34	4	M
SMP-773	Nasal	22/08/2020	10:11	22/08/2020	10:11	+	Nasopharynx	22/08/2020	10:11	22/08/2020	15:00	Bio-tech	+	26,34	4	M
SMP-803	Nasal	24/08/2020	09:55	24/08/2020	09:55	+	Nasopharynx	24/08/2020	09:55	24/08/2020	15:00	Bio-tech	+	26,34	4	M
SMP-704	Nasal	13/08/2020	9:43	13/08/2020	9:43	+	Nasopharynx	13/08/2020	9:43	13/08/2020	15:00	Bio-tech	+	26,36	4	M
SMP-118	Nasal	4/04/2020	09:44	4/04/2020	09:44	+	Nasopharynx	4/04/2020	09:44	4/04/2020	15:00	Bio-tech	+	26,43	5	M
SMP-285	Nasal	4/28/20	08:50	4/28/20	08:50	+	Nasopharynx	4/28/20	08:50	4/28/20	15:00	Bio-tech	+	26,45	4	M
SMP-265	Nasal	4/25/20	08:57	4/25/20	08:57	+	Nasopharynx	4/25/20	08:57	4/25/20	15:00	Bio-tech	+	26,46	2	M
SMP-104	Nasal	3/04/2020	09:03	3/04/2020	09:03	+	Nasopharynx	3/04/2020	09:03	3/04/2020	15:00	Bio-tech	+	26,49	2	M
SMP-48	Nasal	17/03/2020	09:44	17/03/2020	09:44	+	Nasopharynx	17/03/2020	09:44	17/03/2020	15:00	Bio-tech	+	26,66	5	M
SMP-208	Nasal	4/16/20	11:44	4/16/20	11:44	+	Nasopharynx	4/16/20	11:44	4/16/20	15:00	Bio-tech	+	26,66	4	M
SMP-327	Nasal	3/05/2020	11:49	3/05/2020	11:49	+	Nasopharynx	3/05/2020	11:49	3/05/2020	15:00	Bio-tech	+	26,66	3	FM
SMP-360	Nasal	7/05/2020	11:45	7/05/2020	11:45	+	Nasopharynx	7/05/2020	11:45	7/05/2020	15:00	Bio-tech	+	26,66	5	M
SMP-533	Nasal	1/07/2020	10:58	1/07/2020	10:58	+	Nasopharynx	1/07/2020	10:58	1/07/2020	15:00	Bio-tech	+	26,66	4	FM
SMP-605	Nasal	25/07/2020	11:12	25/07/2020	11:12	+	Nasopharynx	25/07/2020	11:12	25/07/2020	15:00	Bio-tech	+	26,66	3	M
SMP-742	Nasal	19/08/2020	11:47	19/08/2020	11:47	+	Nasopharynx	19/08/2020	11:47	19/08/2020	15:00	Bio-tech	+	26,66	3	M
SMP-316	Nasal	1/05/2020	10:59	1/05/2020	10:59	+	Nasopharynx	1/05/2020	10:59	1/05/2020	15:00	Bio-tech	+	26,67	3	M
SMP-511	Nasal	15/06/2020	10:00	15/06/2020	10:00	+	Nasopharynx	15/06/2020	10:00	15/06/2020	15:00	Bio-tech	+	26,74	3	M
SMP-376	Nasal	10/05/2020	8:23	10/05/2020	8:23	+	Nasopharynx	10/05/2020	8:23	10/05/2020	15:00	Bio-tech	+	26,76	3	M
SMP-459	Nasal	5/26/20	11:12	5/26/20	11:12	+	Nasopharynx	5/26/20	11:12	5/26/20	15:00	Bio-tech	+	26,76	3	M
SMP-353	Nasal	7/05/2020	10:11	7/05/2020	10:11	+	Nasopharynx	7/05/2020	10:11	7/05/2020	15:00	Bio-tech	+	26,77	4	M
SMP-760	Nasal	22/08/2020	09:43	22/08/2020	09:43	+	Nasopharynx	22/08/2020	09:43	22/08/2020	15:00	Bio-tech	+	26,77	2	FM
SMP-798	Nasal	24/08/2020	9:50	24/08/2020	9:50	+	Nasopharynx	24/08/2020	9:50	24/08/2020	15:00	Bio-tech	+	26,77	2	M
SMP-516	Nasal	21/06/2020	10:33	21/06/2020	10:33	+	Nasopharynx	21/06/2020	10:33	21/06/2020	15:00	Bio-tech	+	26,78	3	M
SMP-83	Nasal	30/03/2020	10:38	30/03/2020	10:38	+	Nasopharynx	30/03/2020	10:38	30/03/2020	15:00	Bio-tech	+	26,80	2	M
SMP-465	Nasal	5/27/20	11:12	5/27/20	11:12	+	Nasopharynx	5/27/20	11:12	5/27/20	15:00	Bio-tech	+	26,80	3	M
SMP-824	Nasal	26/08/2020	8:46	26/08/2020	8:46	+	Nasopharynx	26/08/2020	8:46	26/08/2020	15:00	Bio-tech	+	26,87	4	M
SMP-84	Nasal	30/03/2020	09:03	30/03/2020	09:03	+	Nasopharynx	30/03/2020	09:03	30/03/2020	15:00	Bio-tech	+	26,88	4	FM
SMP-401	Nasal	5/15/20	09:58	5/15/20	09:58	+	Nasopharynx	5/15/20	09:58	5/15/20	15:00	Bio-tech	+	26,89	4	FM
SMP-106	Nasal	3/04/2020	08:56	3/04/2020	08:56	+	Nasopharynx	3/04/2020	08:56	3/04/2020	15:00	Bio-tech	+	26,98	3	M
SMP-261	Nasal	4/24/20	10:58	4/24/20	10:58	+	Nasopharynx	4/24/20	10:58	4/24/20	15:00	Bio-tech	+	26,98	4	M
SMP-383	Nasal	11/05/2020	09:55	11/05/2020	09:55	+	Nasopharynx	11/05/2020	09:55	11/05/2020	15:00	Bio-tech	+	26,98	4	M
SMP-521	Nasal	24/06/2020	09:21	24/06/2020	09:21	+	Nasopharynx	24/06/2020	09:21	24/06/2020	15:00	Bio-tech	+	26,98	3	M
SMP-326	Nasal	3/05/2020	11:15	3/05/2020	11:15	+	Nasopharynx	3/05/2020	11:15	3/05/2020	15:00	Bio-tech	+	27,01	3	M
SMP-346	Nasal	6/05/2020	09:25	6/05/2020	09:25	+	Nasopharynx	6/05/2020	09:25	6/05/2020	15:00	Bio-tech	+	27,09	3	FM
SMP-306	Nasal	30/04/2020	11:43	30/04/2020	11:43	+	Nasopharynx	30/04/2020	11:43	30/04/2020	15:00	Bio-tech	+	27,11	3	M
SMP-328	Nasal	3/05/2020	12:12	3/05/2020	12:12	+	Nasopharynx	3/05/2020	12:12	3/05/2020	15:00	Bio-tech	+	27,11	5	M
SMP-788	Nasal	24/08/2020	09:02	24/08/2020	09:02	+	Nasopharynx	24/08/2020	09:02	24/08/2020	15:00	Bio-tech	+	27,11	2	M
SMP-834	Nasal	27/08/2020	9:45	27/08/2020	9:45	+	Nasopharynx	27/08/2020	9:45	27/08/2020	15:00	Bio-tech	+	27,11	4	M
SMP-93	Nasal	1/04/2020	10:40	1/04/2020	10:40	+	Nasopharynx	1/04/2020	10:40	1/04/2020	15:00	Bio-tech	+	27,12	2	M
SMP-543	Nasal	6/07/2020	10:44	6/07/2020	10:44	+	Nasopharynx	6/07/2020	10:44	6/07/2020	15:00	Bio-tech	+	27,15	3	M
SMP-312	Nasal	1/05/2020	09:37	1/05/2020	09:37	+	Nasopharynx	1/05/2020	09:37	1/05/2020	15:00	Bio-tech	+	27,18	3	FM
SMP-333	Nasal	4/05/2020	10:44	4/05/2020	10:44	+	Nasopharynx	4/05/2020	10:44	4/05/2020	15:00	Bio-tech	+	27,24	3	M
SMP-613	Nasal	26/07/2020	10:44	26/07/2020	10:44	+	Nasopharynx	26/07/2020	10:44	26/07/2020	15:00	Bio-tech	+	27,27	3	M
SMP-761	Nasal	22/08/2020	09:49	22/08/2020	09:49	+	Nasopharynx	22/08/2020	09:49	22/08/2020	15:00	Bio-tech	+	27,43	3	FM
SMP-835	Nasal	27/08/2020	9:57	27/08/2020	9:57	+	Nasopharynx	27/08/2020	9:57	27/08/2020	15:00	Bio-tech	+	27,44	2	M
SMP-847	Nasal	28/08/2020	13:29	28/08/2020	13:29	+	Nasopharynx	28/08/2020	13:29	28/08/2020	15:00	Bio-tech	+	27,44	3	M
SMP-540	Nasal	5/07/2020	09:24	5/07/2020	09:24	+	Nasopharynx	5/07/2020	09:24	5/07/2020	15:00	Bio-tech	+	27,54	3	M
SMP-260	Nasal	4/24/20	10:24	4/24/20	10:24	+	Nasopharynx	4/24/20	10:24	4/24/20	15:00	Bio-tech	+	27,55	3	M
SMP-10	Nasal	22/03/2020	10:05	22/03/2020	10:05	+	Nasopharynx	22/03/2020	10:05	22/03/2020	15:00	Bio-tech	+	27,66	5	M
SMP-45	Nasal	16/03/2020	11:12	16/03/2020	11:12	+	Nasopharynx	16/03/2020	11:12	16/03/2020	15:00	Bio-tech	+	27,66	3	M
SMP-406	Nasal	5/16/20	10:04	5/16/20	10:04	+	Nasopharynx	5/16/20	10:04	5/16/20	15:00	Bio-tech	+	27,66	4	M
SMP-500	Nasal	5/06/2020	09:22	5/06/2020	09:22	+	Nasopharynx	5/06/2020	09:22	5/06/2020	15:00	Bio-tech	+	27,67	5	M
SMP-710	Nasal	14/08/2020	10:47	14/08/2020	10:47	+	Nasopharynx	14/08/2020	10:47	14/08/2020	15:00	Bio-tech	+	27,70	5	M
SMP-105	Nasal	3/04/2020	08:47	3/04/2020	08:47	+	Nasopharynx	3/04/2020	08:47	3/04/2020	15:00	Bio-tech	+	27,76	4	M
SMP-146	Nasal	8/04/2020	09:03	8/04/2020	09:03	+	Nasopharynx	8/04/2020	09:03	8/04/2020	15:00	Bio-tech	+	27,77	5	FM
SMP-213	Nasal	4/17/20	10:11	4/17/20	10:11	+	Nasopharynx	4/17/20	10:11	4/17/20	15:00	Bio-tech	+	27,77	4	M
SMP-284	Nasal	4/27/20	10:43	4/27/20	10:43	+	Nasopharynx	4/27/20	10:43	4/27/20	15:00	Bio-tech	+	27,77	4	FM
SMP-426	Nasal	5/21/20	09:03	5/21/20	09:03	+	Nasopharynx	5/21/20	09:03	5/21/20	15:00	Bio-tech	+	27,77	5	FM
SMP-488	Nasal	5/30/20	09:44	5/30/20	09:44	+	Nasopharynx	5/30/20	09:44	5/30/20	15:00	Bio-tech	+	27,77	4	M
SMP-566	Nasal	13/07/2020	09:03	13/07/2020	09:03	+	Nasopharynx	13/07/2020	09:03	13/07/2020	15:00	Bio-tech	+	27,77	5	FM
SMP-644	Nasal	30/07/2020	09:03	30/07/2020	09:03	+	Nasopharynx	30/07/2020	09:03	30/07/2020	15:00	Bio-tech				

SMP-46	Nasal	17/03/2020	11:16	17/03/2020	11:16	+	Nasopharynx	17/03/2020	11:16	17/03/2020	15:00	Bio-tech	+	27,88	3	M
SMP-621	Nasal	27/07/2020	09:49	27/07/2020	09:49	+	Nasopharynx	27/07/2020	09:49	27/07/2020	15:00	Bio-tech	+	27,89	3	FM
SMP-75	Nasal	22/03/2020	08:50	22/03/2020	08:50	+	Nasopharynx	22/03/2020	08:50	22/03/2020	15:00	Bio-tech	+	27,91	4	FM
SMP-703	Nasal	13/08/2020	9:13	13/08/2020	9:13	+	Nasopharynx	13/08/2020	9:13	13/08/2020	15:00	Bio-tech	+	27,91	4	FM
SMP-826	Nasal	26/08/2020	9:26	26/08/2020	9:26	+	Nasopharynx	26/08/2020	9:26	26/08/2020	15:00	Bio-tech	+	28,00	4	M
SMP-150	Nasal	8/04/2020	10:42	8/04/2020	10:42	+	Nasopharynx	8/04/2020	10:42	8/04/2020	15:00	Bio-tech	+	28,01	5	M
SMP-827	Nasal	26/08/2020	9:48	26/08/2020	9:48	+	Nasopharynx	26/08/2020	9:48	26/08/2020	15:00	Bio-tech	+	28,01	3	M
SMP-258	Nasal	4/24/20	8:44	4/24/20	8:44	+	Nasopharynx	4/24/20	8:44	4/24/20	15:00	Bio-tech	+	28,09	5	M
SMP-35	Nasal	13/03/2020	08:47	13/03/2020	08:47	+	Nasopharynx	13/03/2020	08:47	13/03/2020	15:00	Bio-tech	+	28,11	4	FM
SMP-580	Nasal	17/07/2020	09:38	17/07/2020	09:38	+	Nasopharynx	17/07/2020	09:38	17/07/2020	15:00	Bio-tech	+	28,14	3	M
SMP-560	Nasal	10/07/2020	09:35	10/07/2020	09:35	+	Nasopharynx	10/07/2020	09:35	10/07/2020	15:00	Bio-tech	+	28,19	3	M
SMP-59	Nasal	21/03/2020	09:25	21/03/2020	09:25	+	Nasopharynx	21/03/2020	09:25	21/03/2020	15:00	Bio-tech	+	28,22	4	FM
SMP-850	Nasal	29/08/2020	9:06	29/08/2020	9:06	+	Nasopharynx	29/08/2020	9:06	29/08/2020	15:00	Bio-tech	+	28,22	3	M
SMP-854	Nasal	29/08/2020	10:32	29/08/2020	10:32	+	Nasopharynx	29/08/2020	10:32	29/08/2020	15:00	Bio-tech	+	28,33	4	M
SMP-750	Nasal	21/08/2020	8:24	21/08/2020	8:24	+	Nasopharynx	21/08/2020	8:24	21/08/2020	15:00	Bio-tech	+	28,45	3	M
SMP-134	Nasal	6/04/2020	10:54	6/04/2020	10:54	+	Nasopharynx	6/04/2020	10:54	6/04/2020	15:00	Bio-tech	+	28,55	4	FM
SMP-577	Nasal	16/07/2020	08:47	16/07/2020	08:47	+	Nasopharynx	16/07/2020	08:47	16/07/2020	15:00	Bio-tech	+	28,61	3	M
SMP-06	Nasal	21/03/2020	10:43	21/03/2020	10:43	+	Nasopharynx	21/03/2020	10:43	21/03/2020	15:00	Bio-tech	+	28,66	5	M
SMP-830	Nasal	26/08/2020	11:00	26/08/2020	11:00	+	Nasopharynx	26/08/2020	11:00	26/08/2020	15:00	Bio-tech	+	28,77	3	M
SMP-70	Nasal	22/03/2020	9:35	22/03/2020	9:35	+	Nasopharynx	22/03/2020	9:35	22/03/2020	15:00	Bio-tech	+	28,88	3	M
SMP-286	Nasal	4/28/20	09:03	4/28/20	09:03	+	Nasopharynx	4/28/20	09:03	4/28/20	15:00	Bio-tech	+	28,88	5	FM
SMP-323	Nasal	2/05/2020	10:58	2/05/2020	10:58	+	Nasopharynx	2/05/2020	10:58	2/05/2020	15:00	Bio-tech	+	28,88	4	FM
SMP-468	Nasal	5/28/20	9:14	5/28/20	9:14	+	Nasopharynx	5/28/20	9:14	5/28/20	15:00	Bio-tech	+	28,88	5	M
SMP-496	Nasal	6/2/20	09:03	6/2/20	09:03	+	Nasopharynx	6/2/20	09:03	6/2/20	15:00	Bio-tech	+	28,88	5	FM
SMP-678	Nasal	7/08/2020	09:44	7/08/2020	09:44	+	Nasopharynx	7/08/2020	09:44	7/08/2020	15:00	Bio-tech	+	28,88	5	M
SMP-36	Nasal	14/03/2020	08:56	14/03/2020	08:56	+	Nasopharynx	14/03/2020	08:56	14/03/2020	15:00	Bio-tech	+	29,00	3	M
SMP-397	Nasal	5/14/20	11:37	5/14/20	11:37	+	Nasopharynx	5/14/20	11:37	5/14/20	15:00	Bio-tech	+	29,00	6	FM
SMP-294	Nasal	4/29/20	09:03	4/29/20	09:03	+	Nasopharynx	4/29/20	09:03	4/29/20	15:00	Bio-tech	+	29,01	6	FM
SMP-856	Nasal	30/08/2020	9:38	30/08/2020	9:38	+	Nasopharynx	30/08/2020	9:38	30/08/2020	15:00	Bio-tech	+	29,01	6	FM
SMP-392	Nasal	13/05/2020	11:25	13/05/2020	11:25	+	Nasopharynx	13/05/2020	11:25	13/05/2020	15:00	Bio-tech	+	29,15	6	M
SMP-587	Nasal	19/07/2020	09:33	19/07/2020	09:33	+	Nasopharynx	19/07/2020	09:33	19/07/2020	15:00	Bio-tech	+	29,29	6	M
SMP-638	Nasal	29/07/2020	9:21	29/07/2020	9:21	+	Nasopharynx	29/07/2020	9:21	29/07/2020	15:00	Bio-tech	+	29,34	7	M
SMP-68	Nasal	22/03/2020	09:44	22/03/2020	09:44	+	Nasopharynx	22/03/2020	09:44	22/03/2020	15:00	Bio-tech	+	29,59	7	M
SMP-142	Nasal	7/04/2020	11:00	7/04/2020	11:00	-	Nasopharynx	7/04/2020	11:00	7/04/2020	15:00	Bio-tech	+	25,70	1	FM
SMP-15	Nasal	23/03/2020	10:38	23/03/2020	10:38	-	Nasopharynx	23/03/2020	10:38	23/03/2020	15:00	Bio-tech	+	27,06	2	M
SMP-389	Nasal	13/05/2020	09:12	13/05/2020	09:12	-	Nasopharynx	13/05/2020	09:12	13/05/2020	15:00	Bio-tech	+	27,87	3	M
SMP-430	Nasal	5/21/20	11:22	5/21/20	11:22	-	Nasopharynx	5/21/20	11:22	5/21/20	15:00	Bio-tech	+	27,91	5	M
SMP-32	Nasal	29/03/2020	09:37	29/03/2020	09:37	-	Nasopharynx	29/03/2020	09:37	29/03/2020	15:00	Bio-tech	+	28,23	3	FM
SMP-98	Nasal	2/04/2020	09:40	2/04/2020	09:40	-	Nasopharynx	2/04/2020	09:40	2/04/2020	15:00	Bio-tech	+	28,33	2	FM
SMP-263	Nasal	4/24/20	11:44	4/24/20	11:44	-	Nasopharynx	4/24/20	11:44	4/24/20	15:00	Bio-tech	+	28,33	3	M
SMP-594	Nasal	21/07/2020	09:03	21/07/2020	09:03	-	Nasopharynx	21/07/2020	09:03	21/07/2020	15:00	Bio-tech	+	28,44	2	M
SMP-03	Nasal	21/03/2020	09:41	21/03/2020	09:41	-	Nasopharynx	21/03/2020	09:41	21/03/2020	15:00	Bio-tech	+	28,55	4	M
SMP-205	Nasal	4/16/20	10:45	4/16/20	10:45	-	Nasopharynx	4/16/20	10:45	4/16/20	15:00	Bio-tech	+	28,81	2	M
SMP-114	Nasal	4/04/2020	10:43	4/04/2020	10:43	-	Nasopharynx	4/04/2020	10:43	4/04/2020	15:00	Bio-tech	+	29,11	6	FM
SMP-472	Nasal	5/28/20	10:33	5/28/20	10:33	-	Nasopharynx	5/28/20	10:33	5/28/20	15:00	Bio-tech	+	29,11	6	FM
SMP-27	Nasal	28/03/2020	09:33	28/03/2020	09:33	-	Nasopharynx	28/03/2020	09:33	28/03/2020	15:00	Bio-tech	+	29,33	7	M
Ct-Average														24,94		
False 25-30 Ct						13										
Total correct						169							182			

Reference kit:

Name: New Coronavirus 2019-nCoV nucleic acid detection kit (fluorescence PCR method)
Registration number: National registration 20203400057
Manufacturer: Shanghai ZJ Bio-Tech Co., Ltd.

Brand Name of VTM System Used:

AS-SWAB UMF AS-VTM Primary Item Number (GTIN): 8683020782808 (GS1) Reference of VTM System: SWPSYS
Ingredients of the VTM system:
• Sterile dacron swab in a plastic package.
• Viral transport medium (20mM Tris (pH 8.0), 2mM EDTA) 3ml in plastic tube.

Antigen test result with the analyzer or reader value:

Novacheck® Antigen Rapid Test is an instrument-free rapid test method. SARS-CoV-2 Antigen Rapid Test employs immuno-lateral chromatography technology for the qualitative detection of antigens. When the antigen advancing in the membrane and reaches the antibody area, it provides a visible color change. Therefore, there is no need to use a device. As a result of the work of the test, the end of the test is visibly interpreted without the need for a device.

After the collection of all the test materials, Rt-PCR tests were performed at the same time at 15:00 every day.